

# Chlorine Disinfection

## Disposable wipes with active chlorine



Item No: 25074



Item No: 25058

**Broad-spectrum disinfection**  
**Easy-to-use**  
**Always correct dosing**

**Wet Wipe Chlorine Disinfection Wipes are impregnated with a granulated chlorine product**

By using Wet Wipe Chlorine Disinfection you will experience a simple work procedure and achieve an optimal disinfection.

The packet contains a burst bag with a surfactant solution. When the bag is burst, the liquid runs in to the wipes, which dissolves the granules that is turned into a solution of active chlorine. The wipes are ready to use after 7 minutes.

**Anwendung**

Use for disinfection of environmental surfaces and contact points on various devices. Disinfects effectively against all viruses for example, Norovirus, HIV and Hepatitis B, vegetative bacteria, and yeast within 1 minute. Contact time against bacterial spores 10 minutes. Can be used for disinfection of surfaces with mold with a contact time of 15 minutes.

Specially recommended for disinfection of contact points for patients in isolation.

The disposable wipes are produced in two sizes; maxi and mini.

**Item no.: 25074: Chlorine Disinfection Maxi**  
Carton of 24 packets  
1 packet contains 5 wipes, 43×30 cm

**Item no.: 25058: Chlorine Disinfection Mini**  
Carton of 20 packets  
1 packet contains 10 wipes, 20×30 cm

**Content in granulated chlorine product**  
Sodium dichloroisocyanurate dihydrate 100 %

**Content in the activated ready to use wipes**  
Active chlorine > 1000–1200 ppm  
Surfactant < 0,5%

**Note**

Activated ready to use wipes retain their disinfection properties up to 7 days after activation. Produced in Germany. Safety data sheet available on request.

**Always read the directions before use.**

Wet Wipe A/S  
info@wetwipe.eu  
Tel: +45 70 266 244  
www.wetwipe.eu

CVR/VAT no: DK 30207726

Wet Wipe GmbH  
info@wetwipe.de  
Tel: +49800 6646940  
www.wetwipe.de

Ust-IdNr. DE 178836011

74-066, EN Ver. 1  
Oct 2017

**wetwipe®**  
Developed in Denmark  
Made in Germany

# Documentation for disinfection effect

## Chlorine Disinfection

Wet Wipe Chlorine Disinfection consists of a structured cloth impregnated with granulated sodium dichloroisocyanurate dihydrate and a burst bag with a surfactant solution. The product is activated by hitting the burst bag, such that the surfactant solution runs into the cloth, which dissolves the granulates that is turned into a solution of hypochlorous acid with an active chlorine concentration of at least 1000 ppm.

Wet Wipe has tested the efficacy of disinfectant solutions of the above granulate under clean conditions with positive results at 1000 ppm active available chlorine. Test results listed in the table below.

**Conclusion:** Incubation time against Clostridium bacterial spores 10 min. Incubation time against vegetative bacteria, yeast fungus, enveloped virus and Norovirus 1 min.

Minimum incubation time for pass of test

Microorganism	EN test	Time
Bacterial spores: <i>Clostridium difficile</i>	EN 13704	10 min
Bacteria: <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i> , <i>Escherichia coli</i> , <i>Pseudomonas aeruginosa</i>	EN 13697	5 min
	EN 1276	1 min
	EN 13727	1 min
	prEN 16615	5 min
Yeast fungus: <i>Candida albicans</i>	EN 13697	5 min
	EN 1650	1 min
	EN 13624	1 min
Mould fungus spores: <i>Aspergillus brasiliensis</i>	EN 13697	15 min
	EN 1650	15 min
	EN 13624	15 min
Virus: <i>Norovirus</i> <i>Adenovirus</i> , <i>poliovirus</i>	EN 14476	1 min
	EN 14476	5 min