

# Ethanol Disinfection

## Disposable wipes with ethanol



Item no.: 51153



Item no.: 51133

**Residue free**  
**Fast acting**  
**Effective against Norovirus**  
**Reclosable click lid**

### Use

Disposable disinfection wipes for disinfection of environmental surface and contact points on devices. Tested according to European test standards.

### Expert reviews

Danish infection hygiene experts at CEI has reviewed Wet Wipe Ethanol Disinfection and concluded that it can be used for disinfection against vegetative bacteria and fungus.

### Storage conditions

Wet Wipes disposable wipes impregnated with ethanol are UN classified as solids containing flammable liquid and may therefore be exempted the technical regulations for storage and handling of flammable liquids.

The disposable wipes are produced in two sizes; maxi and mini.

**Item no.: 51153: Ethanol Disinfection Maxi**  
**Carton of 12 packets**  
**1 packet contains 20 wipes, 43×30 cm**

**Item no.: 51133: Ethanol Disinfection Mini**  
**Carton of 20 packets**  
**1 packet contains 25 wipes, 20×30 cm**

### Content

Ethanol 70,0 %

The ethanol is denaturated with isopropanol (<5%).

Produced in Germany. Safety data sheet available on request.

Always read the directions before use.

Wet Wipe A/S  
info@wetwipe.eu  
Tel: +45 70 266 244  
www.wetwipe.eu

CVR/VAT no: DK 30207726

Wet Wipe GmbH  
info@wetwipe.de  
Tel: +49800 6646940  
www.wetwipe.de

Ust-IdNr. DE 178836011

51-051, EN, Ver. 1  
Oct 2017

**wetwipe®**  
Developed in Denmark  
Made in Germany

# Documentation for disinfection effect

## Ethanol Disinfection

Minimum incubation time for pass of test

Microorganisms	EN test	Exposure time
Bacteria	EN 13727, EN 16615	30 sec
Yeast	EN 13624	30 sec
Mycobacteria	EN 14348	60 sec
Fungus	EN 13624	60 sec
Norovirus	EN 14476	30 sec