

Unides® AF



Ready to use rapid disinfectant based on alcohols for residue free disinfection. Broad antimicrobial spectrum and short drying time. Fragrance free.

Use areas

Use for disinfection of wipeable surfaces and non-invasive alcohol resistant medical devices.

Application

Apply to surface undiluted and allow disinfectant solution to take effect. Keep surface wet for the whole exposure time. Suitable for use with refillable dispenser systems.

Composition

100 g active ingredient solution contains: 24 g Ethanol and 37 g Propan-1-ol.

Antimicrobial Spectrum

Unides® AF is:

- Bactericidal (incl. MRSA)
- Mycobactericidal (*M. terrae*; *M. avium*)
- Yeasticidal, fungicidal (*C. albicans*; *A. brasiliensis*)
- Virus inactivating (HBV/HIV, HCV, BVDV, Rota-, Adeno-, Noro-, Vaccinia-, Influenza- and Herpes viruses)

Listings/registrations

VAH List of Disinfectants

IHO listing

Certified as Medical Device class IIa according to Medical Device Directive 93/42/EEC

BAuA registered according to EU Biocidal Regulation and ECHA guidelines.

Chemical-physical data

Color: colourless
pH-value: 9-10
Flash point: 25°C DIN 51755

Special advice

Use disinfectants safely. Always read the label and product information before use.

4- Field Test

Application in wipe dispenser systems tested in the 4-field test according to EN 16615 and VAH guidelines by Hygieneinstitut IKI - PD Dr.Pitten, Gießen
BAuA Reg.Nr.: N-72471



BAuA Reg.
Nr.: N-72471

Wet Wipe A/S
info@wetwipe.eu
Tel: +45 70 266 244
www.wetwipe.eu

Wet Wipe GmbH
info@wetwipe.de
Tel: +49800 6646940
www.wetwipe.de

52-051, EN Ver. 2
Feb 2018

wetwipe®
Developed in Denmark
Made in Germany

CVR/VAT no: DK 30207726

Ust-IdNr. DE 178836011

Exposure times

	EN test	Exposure times
<i>Bacteria</i>	EN 13727, EN 14561, EN 16615, DGHM/VAH	30 sec
<i>S. aureus (MRSA)</i>	DGHM/VAH	30 sec
<i>Yeast</i>	EN 13624 DGHM/VAH	30 sec
<i>Mycobacteria</i>	EN 14348	60 sec
<i>Fungus</i>	EN 13624 DGHM/VAH	60 sec
<i>HIV/HBV,HCV,Corona</i>	EN 14476	30 sec
<i>Influenza / Herpes</i>	EN 14476	30 sec
<i>BVDV</i>	EN 14476	30 sec
<i>Vaccinia</i>	EN 14476	15 sec
<i>Rotavirus</i>	EN 14476	60 sec
<i>Norovirus</i>	EN 14476	60 sec

Ordering information

	Art.-Nr.	Quantity	Packaging	PZN
Unides AF Liquid				
Spray bottle	50702	250 ml	Carton of 12 bottles	13654708
Unides AF Liquid				
Spray bottle	50204	1.000 ml	Carton of 10 bottles	13654714
Unides AF Liquid				
Canister	50205	5 l		13654720
Canister	50206	10 l		13654737