

# Ethanol Disinfection

## Disposable wipes with ethanol and detergent

Item no.: 25067

Item no.: 25069

**Cleans and disinfects in one step**  
**Fast acting**  
**Effective against Norovirus**  
**Easy-to-use**

### Use

Disposable disinfection wipes for cleaning and disinfection in one step. Use for disinfection of environmental surface and contact points on various devices.

1 min effective incubation time against vegetative bacteria, fungi and virus according to European test standards.

Test report and documentation of efficacy against Norovirus at 30 seconds incubation time is available on request.

### Expert reviews

German and Danish infection hygiene experts at VAH and CEI has reviewed Wet Wipe Ethanol Disinfection and concluded that it can be used for disinfection against bacteria and yeast fungus as well as for disinfection against Norovirus. Expert review available on request.

### Note

Wet Wipes disposable wipes impregnated with ethanol are classified UN 3175 i.e. solids containing flammable liquid and may therefore be exempted the technical regulations for storage and handling of flammable liquids.

The disposable wipes are produced in two sizes; maxi and mini.

Item no.: 25067: Ethanol Disinfection Maxi  
carton of 12 packets  
1 packet contains 20 wipes, 43×30 cm

Item no.: 25069: Ethanol Disinfection Mini  
carton of 20 packets  
1 packet contains 25 wipes, 20×30 cm

### Content

Ethanol 70,0 %  
Detergent < 0,1 %

The ethanol is denaturated with isopropanol. Net content in a package is less than 1 kg, which means the product can be transported without dangerous goods limitation. Carton is therefore marked with LQ label.

Produced in Germany. Safety data sheet available on request.

**Always read the directions before use.**